



New
available with
super coarse
grit.

Time for progress ...

HybridCap

Patent: DE 10 2014 010 033

The (r)evolutionary and modern grit
- **diamond plus ceramic** - combined
with a durable stainless steel cap.

Advantages:

- **Economic efficiency...**
one cap instead of hundreds.
- **Sustainability...**
instead of throwaway.
- **Patient friendliness...**
smoother surfaces in less time.
- **Patient safety and hygiene...**
now you can perform a hygienic and
perfect cleaning, thanks to proper
disinfection or sterilization of your
grinding cap and cap support.




Push[®]

There is no substitute for quality


CE 0124

super coarse hybrid grit • removal of severe callus




Cap and mandrel have to be disinfected and sterilized **unmounted** for hygienic and technical reasons.

coarse grit • removal of thick callus



medium grit • removal of hard callus



Cap support

• long-lasting and exact fitting

fine grit • removal of incipient callus



BUSCH & CO. GmbH
& Co. KG
Unterkaltenbach 17-27
51766 Engelskirchen
GERMANY

Postfach 1152
51751 Engelskirchen
GERMANY

www.busch.eu
mail@busch.eu
Tel.: +49 2263 86-0
Fax: +49 2263 20741



BFLA1/032019